

# Developing a High-Throughput Host Cell Protein Analysis Workflow with Novel SpotMap MS Software


Christopher Sauer<sup>1</sup>, Elsa Gorre<sup>1</sup>, Abby Chiang<sup>2</sup>, Andy Mahan<sup>1</sup>, Will Dracup<sup>3</sup>

(1) Johnson & Johnson Innovative Medicine, Spring House, PA (2) Protein Metrics LLC, Cupertino, CA (3) TotalLab Ltd, Newcastle, United Kingdom


## Introduction

- Host cell proteins (HCPs) are a major concern for biotherapeutic development with impacts to drug efficacy, product shelf life, and adverse reactions.
- The identification and quantitative analysis of host cell proteins is made difficult by the large dynamic range of analytes, especially compared to the biotherapeutic.
- Here we use high speed, narrow window data-independent acquisition (DIA) on a Thermo Orbitrap Astral mass spectrometer to profile HCPs with greater sensitivity, higher throughput, and reduced co-fragmentation of interfering precursor ions.
- By applying HCP analyses earlier in the drug development pipeline, we aim to better understand HCP persistence/clearance and to mitigate risks for downstream purification.

## Orbitrap Astral HCP Optimization



**Figure 2:** Number of unfiltered identifications from different isolation windows. Increasing isolation window width allows more scans and PSMs, but those PSMs are less confident and require more validation.




**Figure 3:** MS/MS data from different isolation windows. Increasing isolation window allows more scans, but loses confidence in IDs due to co-isolation.

**Table 1:** Finalized Parameters for Astral MS


RF Lens (%)	40
AGC Target (%)	500
MS1 Detector	Orbitrap (240k)
MS1 Scan Range	350-1600
Precursor Mass Range	350-1200
MS2 Detector	Astral
MS2 Scan Range	100-1600
HCD Collision Energy (%)	25
Isolation Window	2 Th

- Host cell proteins were analyzed in NIST 8671 reference material.
- Decreasing the scan range and isolation window had the greatest impact on the number of useful spectra collected.
- Methods were optimized to improve **high-quality** identifications.

## Implementation of Routine HCP Analysis



### Longitudinal Study of HCP Profiles



- High risk proteins are flagged at early stages of purification to mitigate risks.
- Tracking of HCP profiles across modalities helps improve understanding of hitchhiking mechanisms

TrendLab used to track individual HCP intensities across multiple samples

## Conclusions

- Native digest followed by LC-MS analysis using high-resolution Orbitrap Astral allows for highly sensitive HCP analyses.
- Methods were optimized to improve the number of useful spectra collected and subsequent HCPs identified.
- SpotMap MS is used for data processing as it is fast and gives the user more control in evaluating potential peptide identifications.
- Future goals are to scale down to nano-flow to improve sensitivity and generate HCP profiles during all cell selection campaigns.

## References

- Huang L, Wang N, Mitchell CE, Brownlee TJ, Maple SR, De MR. A Novel Sample Preparation for Shotgun Proteomics Characterization of HCPs in Antibodies. 2017;89(10):5436-5444. doi:https://doi.org/10.1021/acs.analchem.7b00304
- Damoc E, Tabiwang NA, Pashkova A, Denisov E, Scheffler K, Srzentic K. Increasing sensitivity of high throughput host cell protein analysis on a novel high-resolution accurate mass platform.
- Kiyonami R, Melani R, Chen Y, Leon AD, Du M. Applying UHPLC-HRAM MS/MS Method to Assess Host Cell Protein Clearance during the Purification Process Development of Therapeutic mAbs. International Journal of Molecular Sciences. 2024; 25(17):9687. https://doi.org/10.3390/ijms25179687


## Cell Engineering and Analytical Sciences




Cell Engineering and Analytical Sciences (CEAS) sits at the entry point for biotherapeutics entering the manufacturing pipeline from lead engineering.

- Diverse portfolio of mAbs, multi-specifics, and ADCs

## Current Sample Workflow



**Figure 1:** Native digest is used to minimize interference from the biotherapeutic while still digesting the HCPs



Mobile Phase A: 0.1% FA in H<sub>2</sub>O  
Mobile Phase B: 0.1% FA in ACN  
Flow Rate: 0.3 mL/min




Analytical flow RPLC coupled to Thermo Orbitrap Astral Mass Spectrometer



Proteome Discoverer with Chimerys for improved DIA searching




## SpotMap MS Software



Accession	Priority	Name	Organism	Type	Impact	Trypsin Present	alphaQuant
P55564	<input checked="" type="checkbox"/>	Fructose-bisphosphate aldolase A	Mus musculus			8	
Q91Y99	<input checked="" type="checkbox"/>	Prostaglandin reductase 1	Mus musculus			8	
P40142	<input checked="" type="checkbox"/>	Transferrin	Mus musculus			8	
P02769	<input checked="" type="checkbox"/>	Albumin	Bos taurus			9	

### SpotMap MS from TotalLab

- Quickly visualize TIC from multiple samples to ensure data reproducibility.
- Automatically flag high risk proteins based on published data.



- Designed to rapidly validate identified peptides
- Manual Pass/Fail system based on XIC quality, mass error, and MS/MS fidelity.